

DETERGENTS & CLEANER NU®



NU® GLASS WASHING TABLETS

Perfect hygiene! 500g can with 100 glass detergent tablets for 100 applications. Guaranteed hygienically pure, beer foam-friendly glasses. Already suitably dosed, no quarter of the tablet necessary. Optimum and dry storage in the can.



art-no. 1700410

NU® BRUSH SANITIZER

Disinfection and care! 750g can of brush sanitizer for about 25 applications. Spülboy brush sanitizer ensures constantly fresh brushes and 100% hygiene.



art-no. 1700412

Easy to use!

Add 1 NU® glass detergent tablet into the pre-rinsing canister (brush canister) and, depending on the quantity of glass washes, rinse 100% hygienically for 6-8 hours.

Hint: If no foam be produced, the tablet has completely dissolved. Time to refill another 1 tablet.

Clear advantages!

- Particularly slow release of the NU® tablets guarantees 6-8 hours of **hygienic cleaning power**
- Optimum activation and distribution of hygiene agents through the unique **circulation of fresh water**
- Gentle ingredients of the NU® tablets - **longer life** of the brushes
- **Low foaming** and only! during the active rinsing process. So you have a clear.

Better than liquid!

- **Optimal effect!** NU® tablets do not float on the water, making them last longer.
- **Sustainable hygiene!** NU® glass washing tablets are biodegradable and extremely economical.

Daily Hygiene Plus! NU® Brush cleaner

- To care for and disinfect each Spülboy glass cleaning brush: Carefully remove the brushes, cover with approx. 40 ° C warm water, add 2 parts of clean brush and leave for around 8 hours (overnight). Rinse well in the morning and you're brushes are clean and in a good shape again.